

# Medical Consumables & Personal Protective Equipment

Catalogue



# Contents

<b>Face Mask</b>	<b>Medical Mask</b> 1 - 2 <b>Surgical Mask</b> 3 - 4 <b>Reusable Mask</b> 5 - 6 <b>Particulate Respirator</b> 7 - 20 <b>Protective Mask</b> 21 - 22 <b>Mask Strap</b> 23
<b>Goggles and Face Shields</b>	24 - 26
<b>Gloves</b>	27 - 28
<b>Gowns</b>	29 - 30





EN14683: 2019 TYPE IIR

# Medical Face Mask Single-Use

## PACKAGING SPECIFICATIONS


Units per Box 40 pcs  
Box Dimensions 206 X 103 x 91mm  
Country of Origin PR China



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Outer Material Layer	Spunbond, nonwoven
Middle Filtration Layer	Melt-blown, polypropylene, nonwoven Filtration
Inner Material Layer	Spunbond, nonwoven
Product Dimensions	175 X 95mm
Certifications	EN14683: 2019 TYPE IIR
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Non Sterile



## INSTRUCTIONS FOR USE



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.



Pinch the metal strip to fit the shape of the nose.



Pull bottom of mask below chin.



ASTM F2100-19 Test Report Available

# Medical Mask Non Surgical

## PACKAGING SPECIFICATIONS


Units per Box	50 pcs
Box Dimensions	200 x 100 x 105mm
Boxes per Carton	40 boxes
Carton Dimensions	520 X 420 X 450mm
Carton Gross Weight	9.5kg
Carton Net Weight	5.5kg
Units per Carton	2000 pcs
Country of Origin	PR China



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Outer Material Layer	Spunbond, nonwoven
Middle Filtration Layer	Melt-blown, polypropylene, nonwoven Filtration
Inner Material Layer	Spunbond, nonwoven
Product Dimensions	175 X 95mm
Certifications	ASTM F2100-19 Test Report Available
Storage	Store in a cool, dry, ventilated place
Shelf Life	12 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.





# Surgical Mask Type IIR

EN14683: 2019 TYPE IIR

## PACKAGING SPECIFICATIONS

Units per Box	50 pcs
Box Dimensions	200 x 100 x 90mm
Boxes per Carton	40 boxes
Carton Dimensions	520 X 420 X 390mm
Carton Gross Weight	9.5kg
Carton Net Weight	5.5kg
Units per Carton	2000 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE

EN14683: 2019 TYPE IIR



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part number	PCH- KZM001-001
Outer Material Layer	Melt-blown, polypropylene,
Middle Filtration Layer	nonwoven filtration
	Spunbond, nonwoven
Inner Material Layer	175 X 95mm
Product Dimensions	EN14683: 2019 TYPE IIR
Certifications	Store in a cool, dry,
Storage	ventilated place
Shelf Life	12 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2




Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Surgical Mask Single-Use

EN14683: 2019 TYPE IIR

## PACKAGING SPECIFICATIONS


Units per Box 10 pcs  
Box Dimensions 120 X 125 x 20mm  
Country of Origin PR China



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Outer Material Layer	Spunbond, nonwoven
Middle Filtration Layer	Melt-blown, polypropylene, nonwoven Filtration
Inner Material Layer	Spunbond, nonwoven
Product Dimensions	175 X 95mm
Certifications	EN14683: 2019 TYPE IIR
Storage	Store in a cool, dry, ventilated place
Shelf Life	24 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Antibacterial Washable And Reusable Mask

Sliver ion coated fabric  
Washable and Reusable up to 50 times  
Non Medical. Water Repellent

## COLOR OPTION





## PRODUCT PERFORMANCE

Antibacterial Activity: 99 %  
Effectively Restricts: Staphylococcus Aureus  
Klebsiella Pneumonias  
Escherichia Coli

## PRODUCT SPECIFICATION

Outer Material Layer	Sliver ion coated fabric Antibacterial Water repellent - Repels moisture and droplets.
Middle Filtration Layer	Non-woven fabric
Inner Material Layer	100 % Cotton
Part Number	PCH-MSV-001
Product Dimensions	XS, S, M, L
Certifications	CE, GB32610, FDA, BAT
Storage	Store in a cool, dry, ventilated place

## SIZE OPTION



## INSTRUCTIONS FOR USE



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.



Pull bottom of mask below chin.





# Antibacterial Washable And Reusable Mask

Sliver ion coated fabric  
Washable and Reusable up to 100 times  
Non Medical. Water Repellent

## COLOR OPTION





## PRODUCT PERFORMANCE

Antibacterial Activity: 99 %  
Effectively Restricts: Staphylococcus Aureus  
Klebsiella Pneumonias  
Escherichia Coli

## PRODUCT SPECIFICATION

Outer Material Layer	Sliver ion coated fabric Antibacterial Water repellent - Repels moisture and droplets.
Middle Filtration Layer	Non-woven fabric
Inner Material Layer	90 % Cotton, 10 % Chitosan
Part Number	PCH-MCT-001
Product Dimensions	XS, S, M, L
Certifications	CE, GB32610, FDA, BAT
Storage	Store in a cool, dry, ventilated place

## SIZE OPTION



## INSTRUCTIONS FOR USE



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.



Pull bottom of mask below chin.





# N95 Particulate Respirator

DTC3X

Cool . Comfortable . Collapse Resistant

Headband design style.  
Ultrasonic welding technology, glue-free and odorless


## PACKAGING SPECIFICATION

Units per Box	20 pcs
Boxes per Carton	20 boxes
Carton Dimensions	620 X 290 X 340mm
Carton Gross Weight	6kg
Carton Net Weight	5kg
Country of Origin	PR China



## PRODUCT PERFORMANCE

NIOSH approved: N95  
EN:149 Certified



A pinch-free nose bridge that comfortably contours to the bridge of the nose, eliminating the discomfort and constant re-adjustment that is necessary with metal nosebands. Exclusive Dura-Mesh® shell protects the filter Dura-Meshmedia so it stays cleaner looking longer.  
It resists collapsing in regular use and even under hot and humid conditions.

## Instructions

Failure to follow all instructions and limitations on the use of this respirator and/or failure to wear this respirator during all times of exposure can reduce respirator effectiveness and may result in sickness or death.

Leave the contaminated area immediately if dizziness, irritation or other distress occurs.  
Store the respirator away from contaminated areas when not in use.  
Dispose of used product in accordance with applicable regulations.

## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.  
As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.



EN 149:2001+A1:2009(FFP2NR)

# Single Use Particulate Respirator GL001A

## PACKAGING SPECIFICATIONS


Units per Box 50 pcs  
Box Dimensions 160 X 140 x 185mm  
Country of Origin PR China



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Outer Material Layer	Spunbond, nonwoven
Middle Filtration Layer	Melt-blown, polypropylene, nonwoven Filtration
Inner Material Layer	Spunbond, nonwoven
Product Dimensions	160 X 110mm
Certifications	EN 149:2001+A1:2009(FFP2NR)
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	General Use



## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.  
As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.



# Particulate Respirator

EN149:2001 + A1:2009  
FFP2

## PACKAGING SPECIFICATIONS

Units per Bag	10pcs
Bag Dimensions	200 x 260mm
Carton Dimensions	592 X 402 X 314mm
Carton Gross Weight	6.18kg
Units per Carton	1000pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE

EN 149:2001+A1: 2009  
FFP2



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	MY3D2-001
Outer Material Layer	Nonwoven, 48g/m <sup>2</sup>
Middle Filtration 2nd Layer	Hot air cotton, 40g/m <sup>2</sup>
Middle Filtration 3rd Layer	PTFE Nano fabric 30g/m <sup>2</sup>
Inner Material Layer	Nonwoven, 25g/m <sup>2</sup>
Product Dimension	304 x 106mm
Certifications	EN149:2001+A1:2009
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.  
As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.





# Particulate Respirator

EN149  
FFP2

## PACKAGING SPECIFICATIONS

Units per Box	10 pcs
Box Dimensions	165 x 130 x 100mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 420 X 670mm
Carton Gross Weight	11.2kg
Carton Net Weight	12.2kg
Units per Carton	800 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP2



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9502FC
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration 1st Layer	PP+ES(Hrough-air Nonwoven, 45g/m <sup>2</sup>
Middle Filtration 2nd Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	160 x 105mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.





# Particulate Respirator

EN149  
FFP2

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	165 x 130 x 100mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 420 X 670mm
Carton Gross Weight	16kg
Carton Net Weight	17kg
Units per Carton	1600 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP2



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9502C
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration 1st Layer	PP+ES(Hrough-air Nonwoven, 45g/m <sup>2</sup>
Middle Filtration 2nd Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	160 x 105mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Particulate Respirator

EN149  
FFP2

## PACKAGING SPECIFICATIONS

Units per Box	10 pcs
Box Dimensions	165 x 130 x 100mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 420 X 670mm
Carton Gross Weight	10.6kg
Carton Net Weight	9.6kg
Units per Carton	800 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP2



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9501FL
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	210 x 80mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Particulate Respirator

EN149  
FFP2

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	165 x 130 x 100mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 420 X 670mm
Carton Gross Weight	13.8kg
Carton Net Weight	12.8kg
Units per Carton	1600 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP2



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9501L
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	210 x 80mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Particulate Respirator

EN149  
FFP3

## PACKAGING SPECIFICATIONS

Units per Box	10 pcs
Box Dimensions	165 x 130 x 100mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 420 X 670mm
Carton Gross Weight	14.4kg
Carton Net Weight	15.4kg
Units per Carton	800 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP3



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9902FC
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration 1st Layer	PP+ES(Hrough-air Nonwoven, 45g/m <sup>2</sup>
Middle Filtration 2nd Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	160 x 105mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.





# Particulate Respirator

EN149  
FFP3

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	165 x 130 x 130mm
Boxes per Carton	80 boxes
Carton Dimensions	680 X 540 X 670mm
Carton Gross Weight	22.4kg
Carton Net Weight	23.4kg
Units per Carton	1600 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149  
FFP3



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	YJX-9902C
Outer Material Layer	PP(Nonwoven) 50g/m <sup>2</sup>
Middle Filtration 1st Layer	PP+ES(Hrough-air Nonwoven, 45g/m <sup>2</sup>
Middle Filtration 2nd Layer	PP(Melt-blown) 40g/m <sup>2</sup>
Inner Material Layer	PP(Nonwoven) 30g/m <sup>2</sup>
Product Dimension	160 x 105mm
Certifications	EN149
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pull bottom of mask below chin.

3



Pinch the metal strip to fit the shape of the nose.



# Particulate Respirator

EN149:2001+A1: 2009  
FFP2 NR

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	180 x 75 x 140mm
Boxes per Carton	40 boxes
Carton Dimensions	580 X 375 X 400mm
Carton Gross Weight	7.3kg
Carton Net Weight	8.4kg
Units per Carton	800 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149:2001+A1: 2009  
FFP2 NR



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	AY1680 (Color: White/Green/Blue)
Outer Material Layer	Spunbond nonwoven
Middle Filtration 1st Layer	Skeleton Cotton
Middle Filtration 2nd Layer	Hot air Cotton polypropylene +polypropylene
Middle Filtration 3rd Layer	Meltblown, nonwoven Filtration
Inner Material Layer	Spunbond nonwoven
Product Dimensions	155 X 105mm
Certifications	EN149:2001+A1: 2009 FFP2
Storage	Store in a cool, dry, ventilated place
Shelf Life	24 months



## INSTRUCTIONS FOR USE

1



Wear mask with coloured side facing out, white side facing in.  
Ensure the mask fully covers nose and mouth.  
Loop mask straps around ears.

2



Pinch the metal strip to fit the shape of the  
nose.

3



Pull bottom of mask below chin.



# Particulate Respirator

EN149: 2001+A1: 2009  
FFP2 NR

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	130 x 130 x 205mm
Boxes per Carton	20 boxes
Carton Dimensions	665 X 275 X 425mm
Carton Gross Weight	5kg
Carton Net Weight	6kg
Units per Carton	400 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149:2001+A1: 2009  
FFP2 NR



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Part Number	L-202
Outer Material Layer	140gsm PE
Middle Filtration Layer	50gsm Filter media(P2 filter)
Inner Material Layer	140gsm PE lining
Product Dimensions	130 X 120 X 50mm
Certifications	EN149: 2001+A1: 2009 FFP2
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.  
As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.



# Particulate Respirator

EN149: 2001+A1: 2009  
FFP3 NR

## PACKAGING SPECIFICATIONS

Units per Box	5 pcs
Box Dimensions	130 x 130 x 205mm
Boxes per Carton	20 boxes
Carton Dimensions	715 X 415 X 365mm
Carton Gross Weight	2.7kg
Carton Net Weight	5.7kg
Units per Carton	100 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149:2001+A1: 2009  
FFP3 NR



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	L-203V
Outer Material Layer	75gsm PP
Middle Filtration Layer	100gsm Filter media(P3 filter)
Inner Material Layer	140gsm PP lining
Product Dimensions	200 X 135 X 60mm
Certifications	EN149: 2001+A1: 2009 FFP2
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.

As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.





# Particulate Respirator

EN149: 2001+A1: 2009  
FFP3 NR

## PACKAGING SPECIFICATIONS

Units per Box	12 pcs
Box Dimensions	200 x 100 x 140mm
Boxes per Carton	20 boxes
Carton Dimensions	520 X 420 X 300mm
Carton Gross Weight	5kg
Carton Net Weight	4kg
Units per Carton	240 pcs
Country of Origin	PR China

## PRODUCT PERFORMANCE


EN 149:2001+A1: 2009  
FFP3 NR



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	JTY B-11F
Outer Material Layer	PP Spunbond nonwoven
Middle Filtration 1st Layer	Spunbond, nonwoven
Middle Filtration 2nd Layer	Melt-blown nonwoven filtration
Middle Filtration 3rd Layer	Spunbond, nonwoven
Inner Material Layer	Spunbond nonwoven
Product Dimensions	210 X 100mm
Certifications	EN149: 2001+A1: 2009 FFP3
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE



Block the mask on the chin and pass  
the elastic bands over the head with  
the thumbs and the index fingers



Release the inferior elastic band and  
place it on the neck



Place the superior elastic band on the  
top bottom of the skull



Form the nose clip around the nose  
by pressing firmly on both sides of  
the nose with the fingers



# Particulate Respirator

EN149: 2001+A1: 2009  
FFP2 NR

## PACKAGING SPECIFICATIONS

Units per Box	30 pcs
Box Dimensions	200 x 102 x 140mm
Boxes per Carton	20 boxes
Carton Dimensions	545 X 440 X 360mm
Carton Gross Weight	7.5kg
Carton Net Weight	6.5kg
Units per Carton	600 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE


EN 149:2001+A1: 2009  
FFP2 NR



Discard if  
Damaged



Disposable



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Maximum  
relative humidity of  
storage conditions



## PRODUCT SPECIFICATION

Part Number	JTY A-11
Outer Material Layer	PP Spunbond nonwoven
Middle Filtration 1st Layer	Hot air cotton
Middle Filtration 2nd Layer	Melt-blown nonwoven filtration
Middle Filtration 3rd Layer	Spunbond, nonwoven
Inner Material Layer	Spunbond nonwoven
Product Dimensions	210 X 100mm
Certifications	EN149: 2001+A1: 2009 FFP2
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE



Block the mask on the chin and pass  
the elastic bands over the head with  
the thumbs and the index fingers



Release the inferior elastic band and  
place it on the neck



Place the superior elastic band on the  
top/bottom of the skull



Form the nose clip around the nose  
by pressing firmly on both sides of  
the nose with the fingers



# Protective Face Mask

EN149: 2001+A1: 2009

## PACKAGING SPECIFICATIONS

Units per Box	50 pcs
Box Dimensions	150 x 150 x 160mm
Boxes per Carton	20 boxes
Carton Dimensions	770 X 320 X 340mm
Carton Gross Weight	8.08kg
Carton Net Weight	6.06kg
Units per Carton	1000 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE

EN 149:2001+A1: 2009



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Product Name	NR 952
Outer Material Layer	Spunbond nonwoven
Middle Filtration Layer	Melt-blown nonwoven
Inner Material Layer	Spunbond nonwoven
Product Dimensions	155 X 105mm
Certifications	EN149: 2001+A1: 2009
Storage	Store in a cool, dry, ventilated place
Shelf Life	60 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE

1



Position the respirator under your chin and the nosepiece on the bridge of your nose

2



While holding the respirator in position, pull the TOP headband over your head and rest it at the crown of the back of your head.

As you continue to hold the respirator in position, pull the BOTTOM headband over your head and rest it around your upper neck, below your ears.

3



Using both bands, mold the malleable nosepiece to fit snugly against the nosebridge and face, below the eyes.



# Protective Face Mask

EN149: 2001+A1: 2009

## PACKAGING SPECIFICATIONS

Units per Box	20 pcs
Box Dimensions	190 x 100 x 110mm
Boxes per Carton	20 boxes
Carton Dimensions	415 X 395 X 365mm
Carton Gross Weight	5.36kg
Carton Net Weight	3.6kg
Units per Carton	400 pcs
Country of Origin	PR China



## PRODUCT PERFORMANCE

EN 149:2001+A1: 2009



Discard if Damaged



Disposable



Storage Temperature range of storage conditions



Non-Sterile products




Maximum relative humidity of storage conditions



## PRODUCT SPECIFICATION

Product Name	NR F-820
Outer Material Layer	Spunbond nonwoven
Middle Filtration 1st Layer	Melt-blown nonwoven
Middle Filtration 2nd Layer	Melt-blown nonwoven
Inner Material Layer	Spunbond nonwoven
Product Dimensions	195 X 90mm
Certifications	EN149: 2001+A1: 2009
Storage	Store in a cool, dry, ventilated place
Shelf Life	60 months
Category	Healthcare Worker



## INSTRUCTIONS FOR USE



Block the mask on the chin and pass the elastic bands over the head with the thumbs and the index fingers



Release the inferior elastic band and place it on the neck




Place the superior elastic band on the top bottom of the skull



Form the nose clip around the nose by pressing firmly on both sides of the nose with the fingers





# Mask Strap

Adjustable mask strap for avoidance of ear strain or irritation

## PACKAGING SPECIFICATION

Units per Box	100 pcs
Boxes per Carton	20 boxes
Carton Dimensions	400 X 280 X 265mm
Carton Gross Weight	8.2kg
Country of Origin	PR China

## PRODUCT PERFORMANCE

Size	175 x 24mm
Material	PP + PE
Colour	Translucency Light blue Pink
Weight	3.9g



LATEX



Hand Wash Only



Re Usable



Flexible




Discard if Damaged

AVOID EAR IRRITATION  
SOFT AND FLEXIBLE  
ADJUSTABLE FIT  
RE-USABLE  
WASHABLE - HANDWASH ONLY




## INSTRUCTIONS FOR USE

1




Attach the hooks on the extender to your mask strap loop.

2



Use another strap loop hooks on your extender to attach to the hook loops on the another side of the extender.

3



Pinch the metal strip to fit the shape of the nose.



# Goggles with Indirect Venting


EN166:2001, ANSI Z87.1

## PACKAGING SPECIFICATIONS

Carton Dimensions	550 X 395 X 410mm
Carton Gross Weight	17.5kg
Units per Carton	200 pcs
Country of Origin	PR China

## PRODUCT PERFORMANCE

180° wide vision  
Fully enclosed  
Prevent infection  
Anti-dust  
Anti-fogging  
Anti-splash



Storage Temperature range of storage conditions



Non-Sterile products




Disposable



Discard if Damaged

## PRODUCT SPECIFICATION

Part Number	PCH-DYHF602-001
Material	PVC+PC
Size	190 x 62 x 90mm
Thickness	1.5mm
Weight	80g
Certifications	EN166:2001, ANSI Z87.1
Category	Healthcare worker



## INSTRUCTIONS FOR USE



Place goggles over eyes.



Use the headband to secure the goggles to the head.



Adjust for exact fitting.




# Goggles with Indirect Venting

GB 14866-2006/GB 32166  
EN166:2001, ANSI Z87.1

## PACKAGING SPECIFICATIONS

Carton Dimensions	550 X 395 X 410mm
Carton Gross Weight	17.5kg
Units per Carton	200pcs
Country of Origin	PR China



Storage Temperature range of storage conditions



Non-Sterile products

## PRODUCT PERFORMANCE

180° wide vision  
Fully enclosed  
Prevent infection  
Anti-dust  
Anti-fogging  
Anti-splash




Disposable



Discard if Damaged

## PRODUCT SPECIFICATION

Part Number	PCH-SKLY1A-001
Material	PVC+PC
Size	190 x 62 x 90mm
Thickness	1.5mm
Weight	80g
Certifications	GB 14866-2006/GB 32166 EN166:2001, ANSI Z87.1
Category	Healthcare worker



## INSTRUCTIONS FOR USE



Place goggles over eyes.



Use the headband to secure the goggles to the head.



Adjust for exact fitting.



# Face Shield


GB 14866-2006/GB 32166  
EN 166:2001, ANSI Z87.1

## PACKAGING SPECIFICATIONS

Carton Dimensions	705X 400 X 365mm
Carton Gross Weight	10kg
Units per Carton	200pcs
Country of Origin	PR China

## PRODUCT PERFORMANCE

Prevents splash  
Prevents droplets  
Anti fog



Storage Temperature  
range of  
storage conditions



Non-Sterile  
products




Disposable



Discard if  
Damaged

## PRODUCT SPECIFICATION

Part number	PCH-SKL-M1-001
Material	PET
Size	330 x 220mm
Thickness	0.2mm
Weight	44g
Certifications	GB 14866-2006/GB 32166 EN166:2001, ANSI Z87.1
Category	Healthcare worker



## INSTRUCTIONS FOR USE

1



Tear the double sided protective film off before use face shield

2



Place face shield over face, then use the headband secure the face shield to the head.

3



Adjust for exact fitting.





# Nitrile Gloves

EN 455

## PACKAGING SPECIFICATIONS

Units per Box	100 pcs
Box Dimensions	230 x 125 x 60mm
Boxes per Carton	10 boxes
Carton Dimensions	315 X 258 X 245mm
Units per Carton	1000 pcs
Country of Origin	PR China

## SIZE CHART

Size	Item Code	G.W	N.W
S	NGPF7001	4.0kg	3.0kg
M	NGPF7002	4.5kg	3.5kg
L	NGPF7003	4.8kg	3.8kg
XL	NGPF7004	5.2kg	4.2kg



## PRODUCT PERFORMANCE

Non-sterile  
100 % Latex free



Storage Temperature range of storage conditions



Non-Sterile products



Discard if Damaged



Disposable



## PRODUCT SPECIFICATION

Product Name  
Certifications  
Storage  
Shelf Life  
Category

NGPF 7002  
EN 455  
Store in a cool, dry, ventilated place  
36 months  
Healthcare Worker



Maximum relative humidity of storage conditions



## INSTRUCTIONS FOR USE

1



Grasp the outside of one glove at the wrist. Do not touch your bare skin.

2



Peel off the second glove by putting your fingers inside the glove at the top of your wrist.

3



Turn the second glove inside out while pulling it away from your body, leaving the first glove inside the second.

4



Dispose of the gloves safely. Do not reuse the gloves.



# Vinyl Gloves

EN 455

## PACKAGING SPECIFICATIONS

Units per Box	100 pcs
Box Dimensions	230 x 125 x 60mm
Boxes per Carton	10 boxes
Carton Dimensions	315 X 258 X 245mm
Units per Carton	1000 pcs
Country of Origin	PR China


## SIZE CHART

Size	Item Code	G.W	N.W
S	VGPF7001	5.5kg	4.5kg
M	VGPF7002	6.0kg	5.0kg
L	VGPF7003	6.5kg	5.5kg
XL	VGPF7004	6.8kg	5.8kg



## PRODUCT PERFORMANCE

Non-sterile  
100 % Latex free



Storage Temperature range of storage conditions



Non-Sterile products



Discard if Damaged



Disposable



## PRODUCT SPECIFICATION

Product Name	VGPF 3003
Certifications	EN 455
Storage	Store in a cool, dry, ventilated place
Shelf Life	36 months
Category	Healthcare Worker



Maximum relative humidity of storage conditions



## INSTRUCTIONS FOR USE

1



Grasp the outside of one glove at the wrist. Do not touch your bare skin.

2



Peel off the second glove by putting your fingers inside the glove at the top of your wrist.

3



Turn the second glove inside out while pulling it away from your body, leaving the first glove inside the second.

4



Dispose of the gloves safely. Do not reuse the gloves.



# Disposable Surgical Gown SY1104

EN13795

## PACKAGING SPECIFICATIONS

Material	43g SMS
Color	Blue
Size	L - 122 x 153cm XL - 127 x 172cm
Packaging	1 pc/pouch, 100 pcs/carton
Carton Dimension	52 x 52 x 45cm

## CAUTION

Examine the bag, if seal is broken do not use  
Single use only. Please dispose of the gown properly  
after use  
Use within valid period  
To be kept in clean and dry place

## PRODUCT PERFORMANCE

EN13795 - Pending



Discard if  
Damaged



Disposable



Latex free



EO Sterile

## PRODUCT SPECIFICATION

Part number: PCH-SY1104-001  
Single use  
Anti-fluid  
Anti-alcohol  
Anti-blood  
Anti-static  
Comfortable, Lightweight and Breathable  
Flame retardant  
Non-toxic to skin

## INSTRUCTIONS FOR USE




Allow the body of the gown to fall open.  
Extend arms into the sleeves.



Fasten the velcro.



Fasten inside rear waist ties.



Fasten outside rear waist ties.





# Reinforced Disposable Surgical Gown GA-0308

## PACKAGING SPECIFICATIONS

Units per Bag	1 pc/sterile bag
Units per Carton	30 pcs/carton
Carton Dimension	40 x 26 x 48cm
Country of Origin	PR China

## CAUTION

Examine the bag, if seal is broken do not use  
Single use only. Please dispose of the gown properly  
after use  
Use within valid period  
To be kept in clean and dry place

## PRODUCT PERFORMANCE

CHINA YY0506-2009



Discard if  
Damaged



Disposable



Latex free



EO Sterile

## PRODUCT SPECIFICATION

Part Number	PCH-GA-0308-001
Material	48gsm SMS
Size	L - 130 x 150cm

## COMPOSITION AND FUNCTION

SMS disposable surgical gown with single sterile packing, ultra sonic  
sewed, user from bacteria and liquid. Suitable for doctors in surgery  
room.

## INSTRUCTIONS FOR USE

1



Allow the body of the gown to fall open.  
Extend arms into the sleeves.

2



Fasten the velcro.

3



Fasten inside rear waist ties.

4



Fasten outside rear waist ties.



**Thank you.**